

For real-time processing, the MIUDELLA REMS Healthcare Provider Knowledge Assessment can be completed online at [www.MIUDELLAREMS.com](http://www.MIUDELLAREMS.com). For more information, please call the MIUDELLA REMS at 1-855-337-0772.

## INSTRUCTIONS

### TO BECOME A CERTIFIED HEALTHCARE PROVIDER IN THE REMS AND DISPENSE MIUDELLA FOR INSERTION:

1. Review the MIUDELLA Prescribing Information including the Instructions for Use.
2. Review the **Insertion Training Program for Healthcare Providers**.
3. Complete and submit this **Healthcare Provider Knowledge Assessment** online or print and submit to the REMS via fax to 1-855-781-2610.
4. Complete all required fields on this form to avoid a delay in the enrollment process. You must answer ALL eight questions of this **Healthcare Provider Knowledge Assessment** correctly to become certified.

### COMPLETE ALL REQUIRED FIELDS BELOW:

| HEALTHCARE PROVIDER INFORMATION (*denotes required field) |             |
|-----------------------------------------------------------|-------------|
| *First Name:                                              | *Last Name: |
| *NPI Number:                                              |             |
| *Phone:                                                   | *Email:     |

### KNOWLEDGE ASSESSMENT:

Answer the following questions by selecting the single best answer. You must answer ALL eight questions correctly to become certified in the REMS.

1. You must confirm MIUDELLA is preloaded in the inserter tip before use.
  - a. True
  - b. False
2. When should you twist and remove the green end cap from MIUDELLA?
  - a. When you open the package
  - b. After setting the flange
  - c. After moving the slider to position 2 and confirming you are ready to release MIUDELLA from the inserter
  - d. None of the above
3. After turning the green end cap and exposing MIUDELLA's white release thread, you must pull the release thread completely free from the inserter.
  - a. True
  - b. False
4. The proper way to release MIUDELLA from the inserter is to:
  - a. Cut the white release thread
  - b. Completely remove the white release thread
  - c. Either A or B is acceptable
  - d. None of the above
5. After removing the MIUDELLA inserter, what should you confirm?
  - a. Retrieval threads are not visible
  - b. Retrieval threads are visible past the internal os
  - c. Retrieval threads are visible past the external os
  - d. None of the above
6. MIUDELLA has blue pre-cut retrieval threads.
  - a. True
  - b. False
7. The risk of complications increases when MIUDELLA is not properly inserted.
  - a. True
  - b. False
8. Which of the following are increased risk of complications with improper insertion of MIUDELLA?
  - a. Perforation or expulsion
  - b. Infection
  - c. Unintended pregnancy or pregnancy loss
  - d. All of the above

**If you have questions about or require assistance with the MIUDELLA REMS,  
call 1-855-337-0772 Monday-Friday, 8:00 AM – 8:00 PM ET.**

**Report adverse events to Organon LLC at 1-844-674-3200 or  
to the FDA at 1-800-FDA-1088 or online at [www.fda.gov/medwatch](http://www.fda.gov/medwatch).**